

TWF2 Rabbit pAb

货号: B23501

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IP
推荐浓度	IP: 50 - 1:200
理论分子量	39kDa
实测分子量	39kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cell projection,Cytoplasm,cytoskeleton,perinuclear region,stereocilium
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-369 of human TWF2 (NP_009215.1).
序列	MAHQTGIHATEELKEFFAKARAGSVRLIKVVIIEDEQLVLGASQEPVGRWDQDYDRAVLPLLDAQQPCYLLYRLDSQNAQG FEWLFLAWSPDNSPVRLKMLYAATRATVKKEFGGGHIKDELFGTVKDDLSFAGYQKHLSSCAAPAPL TSAERELQQIRINE VKTEISVESKHQTLQGLAFPLQPEAQRALQQLKQKMVNYIQMKDLERETIELVHTEPTDVAQLPSRVPRDAARYHFFLYKH THEGDPLESVVFIYSMPGYKCSIKERMLYSSCKSRLLDSVEQDFHLEIAKKIEIGDGAELTAEFLYDEVHPKQHAFKQAF AKP KPGGKRGHKRLIRGPGENGDDS

靶点信息

研究背景	The protein encoded by this gene was identified by its interaction with the catalytic domain of protein kinase C-zeta. The encoded protein contains an actin-binding site and an ATP-binding site. It is most closely related to twinfilin (PTK9), a conserved actin monomer-binding protein.
------	---

基因ID	11344
基因名	TWF2
Swiss	Q6IBS0
别名	TWF2; A6RP; A6r; MSTP011; PTK9L; twinfilin-2

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn