

HERC4 Rabbit pAb

货号: B23534

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	12kDa/13kDa/106kDa/109kDa/117kDa/118kDa
实测分子量	12kDa/13kDa/106kDa/109kDa/117kDa/118kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cytoplasm,cytosol
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 400-700 of human HERC4 (NP_056416.2).
序列	WTVNEALIQKWSYPSGRFPVEIANEIDGTFSSSGCLNGSFLAVSNDHDHYRTGTRFSGVDMNAARLLFHKLIQPDHPQIS QQVAASLEKNLIPKLTSSLPDVEALRFYLTLPCEPLMSDSNNFTTIAIPFGTALVNLEKAPLKVLENWWSVLEPPLFLKIVELFK EVVVHLLKLYKIGIPPSERRIFNSFLHTALKVLEILHRVNEKMGQIIQYDKFYIHEVQELIDIRNDYINWVQQQAYGMLADIPVTI CTYPFVFDAQAKTTLLQTDAVLQMOMAIDQHRQNVSSLFLPVIESVNPC

靶点信息

研究背景	HERC4 belongs to the HERC family of ubiquitin ligases, all of which contain a HECT domain and at least 1 RCC1 (MIM 179710)-like domain (RLD). The 350-amino acid HECT domain is predicted to catalyze the formation of a thioester with ubiquitin before transferring it to a substrate, and the RLD is predicted to act as a guanine nucleotide exchange factor for small G proteins (Hochrainer et al., 2005 [PubMed 15676274]).[supplied by OMIM, Mar 2008]
------	--

基因ID	26091
基因名	HERC4
Swiss	Q5GLZ8
别名	HERC4; probable E3 ubiquitin-protein ligase HERC4

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn