

PFKP Rabbit pAb

货号: B23588

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	85kDa
实测分子量	85kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 605-784 of human PFKP (NP_002618.1).
序列	IFEEPFDIRDLQSNVEHLTEKMKTTIQRGLVLRNESCSENYTTDFIYQLYSEEGKGVFDCRKNVLGHMQQGGAPSPFDRNF GTKISARAMEWITAKLKEARGRGKKFTTDDSIKVLGSKRNVIFQPVAELKKQTDFEHRIPKEQWWLKLRLMKILAKYKAS YDVSDSGQLEHVQPWSV

靶点信息

研究背景	This gene encodes a member of the phosphofructokinase A protein family. The encoded enzyme is the pl atelet-specific isoform of phosphofructokinase and plays a key role in glycolysis regulation. This gene ma y play a role in metabolic reprogramming in some cancers, including clear cell renal cell carcinomas, and cancer of the bladder, breast, and lung. Alternative splicing results in multiple transcript variants.
基因ID	5214

基因名	PFKP
Swiss	Q01813
别名	PFKP; ATP-PFK; PFK-C; PFK-P; PFKF; phosphofructokinase; platelet

产品验证



Western blot analysis of PFKP expressed in CACO2 using PFKP Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn