

USP4 Rabbit pAb

货号: B23633

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	35kDa/103kDa/108kDa
实测分子量	35kDa/103kDa/108kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat cells, NCI-H1299 cells, HEK-293 cells, RAW 264.7 cells
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 560-700 of human USP4 (NP_003354.2).
序列	DDIFVYEV CSTSV DGS ECV TLPVYFRERKSRPSSTSSASALYGQPLLLSVPKHKLTLES LYQAVCDRISRYVKQPLPDEFGS SPLEPGACNGSRNSCEGEDEEEMEHQEEGKEQLSETEGSGEDEPGNDPSETTQKKIKGQ

靶点信息

研究背景	The protein encoded by this gene is a protease that deubiquitinates target proteins such as ADORA2A and TRIM21. The encoded protein shuttles between the nucleus and cytoplasm and is involved in maintaining operational fidelity in the endoplasmic reticulum. Three transcript variants encoding different isoforms have been found for this gene.
基因ID	7375

基因名	USP4
Swiss	Q13107
别名	USP4; UNP; Unph; ubiquitin specific peptidase 4

产品验证



Western blot analysis of USP4 expressed in Jurkat cells,NCI-H1299 cells,HEK-293 cells using USP4 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug p er lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn