

PROM1 Rabbit pAb

货号: B23659

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IP IF/ICC
推荐浓度	WB: 1:100 - 1:500 IHC: 1:50 - 1:200 IP: 1:500 - 1:1000 IF/ICC: 1:50 - 1:200
理论分子量	92kDa/93kDa/94kDa/96kDa/97kDa
实测分子量	120KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,HCT116
细胞定位	Apical cell membrane,Cell projection,Endoplasmic reticulum,Endoplasmic reticulum-Golgi intermediate co mpartment,Multi-pass membrane protein,cilium,microvillus membrane,photoreceptor outer segment
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 179-433 of human PRO M1 (NP_006008.1).
序列	ANHQVRTRIKRSRKLADSNFKDLRTLNETPEQIKYILAQYNTTKDKAFTDLNSINSVLGGGILDRLRPNIIPVLDEIKSMATAI KETKEALENMNSTLKSLHQSTQLSSSLTSVKTSRSSLNDPLCLVHPSSETCNSIRLSLSQLNSNPELRQLPPVDAELDNV NNVLRTDLDGLVQQGYQSLNDIPDRVQRQTTTVVAGIKRVLNSIGSDIDNVQRLPIQDILSAFSVYVNNTESYIHRNLPTL EEYDSY

靶点信息

研究背景	This gene encodes a pentaspan transmembrane glycoprotein. The protein localizes to membrane protrusions and is often expressed on adult stem cells, where it is thought to function in maintaining stem cell properties by suppressing differentiation. Mutations in this gene have been shown to result in retinitis pigmentosa and Stargardt disease. Expression of this gene is also associated with several types of cancer. This gene is expressed from at least five alternative promoters that are expressed in a tissue-dependent manner. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	8842
基因名	PROM1
Swiss	O43490
别名	PROM1; AC133; CD133; CORD12; MCDR2; MSTP061; PROML1; RP41; STGD4; prominin-1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn