

Podoplanin Rabbit pAb

货号: B23685

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200
理论分子量	12kDa/16kDa/18kDa/24kDa
实测分子量	12kDa/16kDa/18kDa/24kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cell projection,Membrane,Single-pass type I membrane protein,filopodium membrane,lamellipodium membrane,microvillus membrane,ruffle membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 99-199 of human Podoplanin (NP_006465.3).
序列	ASTGQPEDDTETTGLEGGVAMPGAEDDVVTPGTSEDRYKSLTTLVATSVNSVTGIRIEDLPTSESTVHAQEQSPSATASNVATSHSTEKVDGDTQTTVEK

靶点信息

研究背景	This gene encodes a type-I integral membrane glycoprotein with diverse distribution in human tissues. The physiological function of this protein may be related to its mucin-type character. The homologous protein in other species has been described as a differentiation antigen and influenza-virus receptor. The specific function of this protein has not been determined but it has been proposed as a marker of lung injury. Alternatively spliced transcript variants encoding different isoforms have been identified.
------	---

基因ID	10630
基因名	PDPN
Swiss	Q86YL7
别名	PDPN; AGGRUS; GP36; GP40; Gp38; HT1A-1; OTS8; PA2.26; T1A; T1A-2; T1A2; TI1A; podoplanin

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn