

IFNAR2 Rabbit pAb

货号: B23692

产品信息

反应	Mouse,Rat
宿主	rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:100 - 1:500
理论分子量	27kDa/37kDa/57kDa
实测分子量	58KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Rat liver
细胞定位	Membrane,Secreted,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 370-475 of human IFNA R2 (NP_997468.1).
序列	LPEVDVELPTMPKDSPQQLELLSGPCERRKSPLQDPFPEEDYSSTEGSGGRITFNVDLNSVFLRVLDDEDSDDLEAPLMLS SHLEEMVDPEDPDNVQSNHLLASGE

靶点信息

研究背景	The protein encoded by this gene is a type I membrane protein that forms one of the two chains of a receptor for interferons alpha and beta. Binding and activation of the receptor stimulates Janus protein kinase s, which in turn phosphorylate several proteins, including STAT1 and STAT2. Multiple transcript variants encoding at least two different isoforms have been found for this gene.
基因ID	3455

基因名	IFNAR2
Swiss	P48551
别名	IFNAR2; IFN-R; IFN-alpha-REC; IFNABR; IFNARB; IMD45; interferon alpha/beta receptor 2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn