

JMJD6 Rabbit pAb

货号: B23737

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IHC
推荐浓度	IHC: 1:50 - 1:200
理论分子量	43kDa/46kDa/47kDa
实测分子量	43kDa/46kDa/47kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Nucleus,nucleolus,nucleoplasm
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-320 of human JMJD6 (NP_055982.2).
序列	MNHKSKKRIREAKRSARPELKDSLWTRHNYYESFSLSPAADVADNVERADALQLSVEEFVERYERPYKPVVLLNAQEGWS AQEKWTLERLKRKYRNQKFKCGEDNDGYSVKMKMKYIEYMESTRDDSPYIFDSSYGEHPKRRKLLDYKVPKFFTDLL FQYAGEKRRPPYRWFVMGPPRSGTGIHIDPLGTSAWNALVQGHKRWCLFPTSTPRELIKVTRDEGGNQQDEAITWFNVIY PRTQLPTWPPEFKPLEILQKPGETVFPVGGWWHVVLNLDTTIAITQNFASSTNFPVVVHKTVRGRPKLSRKWYRILKQEH

靶点信息

研究背景	This gene encodes a nuclear protein with a JmjC domain. JmjC domain-containing proteins are predicted to function as protein hydroxylases or histone demethylases. This protein was first identified as a putative phosphatidylserine receptor involved in phagocytosis of apoptotic cells; however, subsequent studies have indicated that it does not directly function in the clearance of apoptotic cells, and questioned whether it is a true phosphatidylserine receptor. Multiple transcript variants encoding different isoforms have been found for this gene.
------	---

基因ID	23210
基因名	JMJD6
Swiss	Q6NYC1
别名	JMJD6; PSR; PTDSR; PTDSR1; arginine demethylase and lysine hydroxylase

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn