

MICU1 Rabbit pAb

货号: B23782

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:100 - 1:500
理论分子量	
实测分子量	46KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	OVCAR3,A-431,293T
细胞定位	mitochondrial inner membrane,mitochondrial intermembrane space,mitochondrial membrane,mitochondrion
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 105-255 of human MICU 1 (NP_001182447.1).
序列	GFRDRKVMHEYENRIRAYSTPDKIFRYFATLKVISEPGAEVFMTPEDFVRSITPNEKQPEHLGLDQYIIKRFDGKKISQEREKFADEGSIFYTLGECGLISFSDYIFLTTLVSTPQRNFEIAFKMFDLNGDGEVDMEEFEQVQSIIRSQTS

靶点信息

研究背景	This gene encodes an essential regulator of mitochondrial Ca2+ uptake under basal conditions. The encoded protein interacts with the mitochondrial calcium uniporter, a mitochondrial inner membrane Ca2+ channel, and is essential in preventing mitochondrial Ca2+ overload, which can cause excessive production of reactive oxygen species and cell stress. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Mar 2013]
------	---

基因ID	10367
基因名	MICU1
Swiss	Q9BPX6
别名	CALC;CBARA1;EFHA3;MPXPS;MICU1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn