

NG2/CSPG4 Rabbit pAb

货号: B23833

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	250kDa
实测分子量	251KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,MCF7,Mouse heart
细胞定位	apical plasma membrane,cell surface,extracellular exosome,extracellular region,focal adhesion,Golgi lumen,lamellipodium membrane,lysosomal lumen,nucleoplasm,plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1340-1563 of human NG2/CSPG4 (NP_001888.2).
序列	LGAPLEGVLVELEVLPAAIPLAQNFVSPEGGSLTAPLLRVSGPYFPTLLGLSLQVLEPPQHGAHQKEDGPQARTLSAFS WRMVEEQLRIRYVHDGSETLTDSFVLMANASEMDRQSHPVAFVTVLPVNDQPPILTTNTGLQMWEGATAPIPAEALRSTD GDSGSEDLVYTIEQPSNGRVVLRGAPGTEVRSFTQAQLDGGVLVFSHRGTLDGGFRFRLSDG

靶点信息

研究背景	A human melanoma-associated chondroitin sulfate proteoglycan plays a role in stabilizing cell-substratum interactions during early events of melanoma cell spreading on endothelial basement membranes. CSPG4 represents an integral membrane chondroitin sulfate proteoglycan expressed by human malignant melanoma cells.
------	---

基因ID	1464
基因名	CSPG4
Swiss	Q6UVK1
别名	CSPG4;HMW-MAA;MCSP;MCSPG;MEL-CSPG;MSK16;NG2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn