

SFRP2 Rabbit pAb

货号: B23836

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IF/ICC: 1:50 - 1:200
理论分子量	33kDa
实测分子量	33KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	extracellular space
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SFRP2 (NP_003004.1).
序列	MLQGPGSLLLLFLASHCCLGSARGLFLFGQPDFS YKRSNCKPIPANLQLCHGIEYQNMRLPNLLGHETMKEVLEQAGAWIP LVMKQCHPDTKKFLCSLFA

靶点信息

研究背景	This gene encodes a member of the SFRP family that contains a cysteine-rich domain homologous to the putative Wnt-binding site of Frizzled proteins. SFRPs act as soluble modulators of Wnt signaling. Methylation of this gene is a potential marker for the presence of colorectal cancer.
基因ID	6423

基因名	SFRP2
Swiss	Q96HF1
别名	SFRP2;FRP-2;SARP1;SDF-5

产品验证



Western blot analysis of SFRP2 expressed in HeLa using SFRP2 Rabbit pAb at 1:1000. Secondary antibody : HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn