

INTS11 Rabbit pAb

货号: B23854

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:100 - 1:500
理论分子量	56kDa/64kDa/65kDa/67kDa/68kDa
实测分子量	68KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293F
细胞定位	blood microparticle,cytoplasm,cytosol,nucleoplasm,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 331-600 of human INTS 11 (NP_060341.2).
序列	GQSLQIFRKWAGNEKNMVIMPGYCVQGTVGHKILSGQRKLEMEGRQVLEVKMQVEYMSFSAHADAKGIMQLVQGAEPE SVLLVHGEAKKMEFLKQKIEQELRVNCYPANGETVTLPTSPSIPVGISLGLLKREMAQGLLPEAKKPRLLHGTLMKDSNFR LVSSEQALKELGLAEHQLRFTCRVHLHDTRKEQETALRVYSHLKSVLKDHCVQHLPDGSVTVESVLLQAAAPSEDPGTKVL LVSWTYQDEELGSFLTSLKKGLPQAPS

靶点信息

研究背景	The Integrator complex contains at least 12 subunits and associates with the C-terminal domain of RNA polymerase II large subunit (POLR2A; MIM 180660) and mediates the 3-prime end processing of small nuclear RNAs U1 (RNU1; MIM 180680) and U2 (RNU2; MIM 180690). INTS11, or CPSF3L, is the catalytic subunit of the Integrator complex (Baillat et al., 2005 [PubMed 16239144]).[supplied by OMIM, Mar 2008]
------	---

基因ID	54973
基因名	INTS11
Swiss	Q5TA45
别名	INTS11;CPSF3L;CPSF73L;INT11;RC-68;RC68

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn