

Munc 13-1 Rabbit pAb

货号: **B24261**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:1000 - 1:5000
理论分子量	193kDa
实测分子量	220KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	neuromuscular junction,neuron projection,plasma membrane,presynaptic active zone,presynaptic active zone cytoplasmic component,presynaptic membrane,synaptic vesicle membrane,terminal bouton
纯化	Affinity purification

抗原信息

抗原信息	Recombinant protein of human Munc 13-1.
序列	

靶点信息

研究背景	This gene encodes a member of the UNC13 family. UNC13 proteins bind to phorbol esters and diacylglycerol and play important roles in neurotransmitter release at synapses. Single nucleotide polymorphisms in this gene may be associated with sporadic amyotrophic lateral sclerosis.
基因ID	23025
基因名	UNC13A

Swiss	Q9UPW8
别名	Munc13-1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn