

CD47 Rabbit pAb

货号: B24644

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	33kDa
实测分子量	45-60KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse thymus,Rat spleen,Rat lung
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant protein of mouse CD47.
序列	

靶点信息

研究背景	Predicted to enable protein binding activity involved in heterotypic cell-cell adhesion and thrombospondin receptor activity. Involved in regulation of nitric oxide biosynthetic process. Acts upstream of or within several processes, including monocyte aggregation; opsonization; and positive regulation of phagocytosis . Located in extracellular exosome. Is integral component of plasma membrane. Is expressed in several structures, including alimentary system; central nervous system; genitourinary system; hemolymphoid system gland; and liver and biliary system. Orthologous to human CD47 (CD47 molecule).
基因ID	16423

基因名	
Swiss	Q61735
别名	IAP; Itgp; 9130415E20Rik; B430305P08Rik

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn