

# CD45RB Rabbit pAb

货号: B24785

## 产品信息

|       |  |
|-------|--|
| 反应    | Human  |
| 宿主    | Rabbit   |
| 克隆性   | Polyclonal   |
| 预测反应  |  |
| 应用    | WB IF/ICC  |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:1000<br><b>IF/ICC:</b> 1:50 - 1:200   |
| 理论分子量 | 135kDa   |
| 实测分子量 | 200-250KDa   |
| 形式    | Liquid   |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated   |
| 阳性对照  | Jurkat   |
| 细胞定位  | Membrane,Membrane raft,Single-pass type I membrane protein   |
| 纯化    | Affinity purification  |


## 抗原信息

|      |   |
|------|---|
| 抗原信息 | Recombinant protein of human CD45RB.  |
| 序列   | GMVSTFEQYQFLYDVIASTYPAQNGQVKKNNHQEDKIEFDNEVDKVKQDANCVNPLGAPEKLPEAKEQAEGSEPTSGTE<br>GPEHSVNGPASPALNQGS |

## 靶点信息

|       |  |
|-------|--|
| 研究背景  | The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitosis, and oncogenic transformation. This PTP contains an extracellular domain, a single transmembrane segment and two tandem intracytoplasmic catalytic domains, and thus is classified as a receptor type PTP. This PTP has been shown to be an essential regulator of T- and B-cell antigen receptor signaling. It functions through either direct interaction with components of the antigen receptor complexes, or by activating various Src family kinases required for the antigen receptor signaling. This PTP also suppresses JAK kinases, and thus functions as a regulator of cytokine receptor signaling. Alternatively spliced transcripts variants of this gene, which encode distinct isoforms, have been reported. |
| 基因ID  | 5788   |
| 基因名   | PTPRC  |
| Swiss | P08575-9   |
| 别名    | B220;CD45;CD45R;GP180;L-CA;LCA;LY5;T200;PTPRC  |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)