

CD107b Rabbit pAb

货号: B24786

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	45kDa
实测分子量	102KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,NIH/3T3
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant protein of mouse CD107b.
序列	

靶点信息

研究背景	Enables protein domain specific binding activity. Involved in several processes, including autophagosome maturation; protein stabilization; and protein targeting to lysosome involved in chaperone-mediated autophagy. Acts upstream of or within muscle cell cellular homeostasis. Located in late endosome membrane; lysosomal membrane; and phagocytic vesicle membrane. Is integral component of autophagosome membrane. Is expressed in central nervous system; egg cylinder; lung; and retina. Used to study Danon disease. Human ortholog(s) of this gene implicated in Danon disease and hypertrophic cardiomyopathy. Orthologous to human LAMP2 (lysosomal associated membrane protein 2).
基因ID	16784

基因名	
Swiss	P17047
别名	Mac3; LGP-B; CD107b; Lamp-2; Lamp II; Lamp-2a; Lamp-2b; Lamp-2c

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn