

NLRC4 Rabbit pAb

货号: B25047

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IF: HUVEC cells , Rattus norvegicus WB: N9 cells , Mus musculus , Homo sapiens , Rattus norvegicus IHC: Homo sapiens
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IF/ICC: 1:50 - 1:200
理论分子量	10kDa/18kDa/40kDa/116kDa
实测分子量	110KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse large intestine
细胞定位	Cytoplasm,cytosol
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 85-170 of human NLRC4 (NP_001186067.1).
序列	NGQSLFHQTSEGLDDLAQDLKDLYHTPSFLNFYPLGEDIDIIFNLKSTFTEPVLWRKDQHHRVEQLTLNGLLQALQSPC IIEGE

靶点信息

研究背景	This gene encodes a member of the caspase recruitment domain-containing NLR family. Family members play essential roles in innate immune response to a wide range of pathogenic organisms, tissue damage and other cellular stresses. Mutations in this gene result in autoinflammation with infantile enterocolitis. A lternative splicing results in multiple transcript variants.
------	--

基因ID	58484
基因名	NLRC4
Swiss	Q9NPP4
别名	NLRC4;AIFEC;CARD12;CLAN;CLAN1;CLANA;CLANB;CLANC;CLAND;CLR2.1;FCAS4;IPAF

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn