

GSDMD (Full Length+C terminal) Rabbit pAb

货号: B25049

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens , Mus musculus
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IF/ICC: 1:50 - 1:200
理论分子量	53kDa
实测分子量	53KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-431,THP-1
细胞定位	cytosol,extracellular region,extracellular space,NLRP3 inflammasome complex,nucleoplasm,plasma mem brane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 276-484 of human GSD MD (NP_079012.3).
序列	GVPAEGAFTEDFQGLRAEVETISKELELLDRELCQLLLEGLEGVLRDQLALRALEEALQGGSLGPVEPLDGPAGAVLECLV LSSGMLVPELAIPVVYLLGALTMLSETQHKLLAEALESQTLLGPLELVGSLLEQSAPWQERSTMSLPPGLLGNSWGEGAPA WVLLDECGLELGEDTPHVCWEPQAQGRMCALYASLALLSGLSQEPH

靶点信息

研究背景	Gasdermin D is a member of the gasdermin family. Members of this family appear to play a role in regula tion of epithelial proliferation. Gasdermin D has been suggested to act as a tumor suppressor. Alternative ly spliced transcript variants have been described.
------	--

基因ID	79792
基因名	GSDMD
Swiss	P57764
别名	GSDMD;DF5L;DFNA5L;FKSG10;GSDMDC1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn