

Phospho-Syk-Y525/526 Rabbit pAb

货号: B25069

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB IP
推荐浓度	WB: 1:500 - 1:2000 IP: 1:50 - 1:100
理论分子量	69kDa/72kDa
实测分子量	72kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat
细胞定位	Cell membrane,Cytoplasm,cytosol
纯化	Affinity purification

抗原信息


抗原信息	A synthetic phosphorylated peptide around Y525 & Y526 of human Syk (NP_003168.2).
序列	ENYYKA

靶点信息

研究背景	This gene encodes a member of the family of non-receptor type Tyr protein kinases. This protein is widely expressed in hematopoietic cells and is involved in coupling activated immunoreceptors to downstream signaling events that mediate diverse cellular responses, including proliferation, differentiation, and phagocytosis. It is thought to be a modulator of epithelial cell growth and a potential tumour suppressor in human breast carcinomas. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
基因ID	6850

基因名	SYK
Swiss	P43405
别名	SYK;p72-Syk

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn