

Rig-I/DDX58 Rabbit pAb

货号: B25072

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	101kDa/106kDa
实测分子量	102-110kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Mouse brain,Mouse ovary,Rat brain
细胞定位	Cell junction,Cell projection,Cytoplasm,cytoskeleton,ruffle membrane,tight junction
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 726-925 of human Rig-I/DDX58 (NP_055129.2).
序列	QTRGRGRARGSKCFLLTSNAGVIEKEQINMYKEKMMNDSILRLQTWDEAVFREKILHIQTHEKFIRDSQEKPKPVPDKENK KLLCRKCKALACYTADVRVIEECHYTVLGDAFKECFVSRPHPKPKQFSSF EKRAKIFCARQNC SHDWGIHVKYKTFEIPVIKIE SFVVEDIATGVQTLYSKWKDFHFEKIPFDP AEMSK

靶点信息

研究背景	DEAD box proteins, characterized by the conserved motif Asp-Glu-Ala-Asp (DEAD), are putative RNA helicases which are implicated in a number of cellular processes involving RNA binding and alteration of RNA secondary structure. This gene encodes a protein containing RNA helicase-DEAD box protein motifs and a caspase recruitment domain (CARD). It is involved in viral double-stranded (ds) RNA recognition and the regulation of immune response.
------	---

基因ID	23586
基因名	DDX58
Swiss	O95786
别名	DDX58;RIG-I;RIGI;RLR-1;SGMRT2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn