

β-Actin Rabbit pAb

货号: B25118

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: MSC cells , Mus musculus , Rattus norvegicus , Arabidopsis thaliana , Homo sapiens , Tree , Flowers , Melodinus cochinchinensis , Sus scrofa , flounder , fish , lamprey , Gallus gallus , Crassostrea gigas , acy g IF: Homo sapiens,Rattus norvegicus , Homo sapiens , Mus musculus IP: Mus musculus , Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	41kDa
实测分子量	42kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,MCF7,A-431,NIH/3T3,C6,Mouse testis,Mouse kidney,Rat testis
细胞定位	Cytoplasm,cytoskeleton
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-375 of human β-Actin (NP_001092.1).
序列	MDDDDIAALVVDNGSGMCKAGFAGDDAPRAVFPSIVGRPRHQGVMVGMGQKDSYVGDEAQSKRGILTLKYPIEHGIVTN WDDMEKIWHHTFYNELRVAPEEHPVLLTEAPLNPKANREKMTQIMFETFNTPAMYVAIQAVLSLYASGRTTGIVMDSGDG VTHTVPIYEGYALPHAILRLDLAGRDLDYLMKILTERGYSFTTTAEREIVRDIKEKLCYVALDFEQEMATAASSSSLEKSYELP DGQVITIGNERFRCPEALFQPSFLGMESCGIHETTFNSIMKCDVDIRKDLYANTVLSGGTTMYPGIADRMQKEITALAPSTM KIKIIPPERKYSVWIGGSILASLSTFQQMWISKQEYDESGPSIVHRKCF

靶点信息

研究背景	This gene encodes one of six different actin proteins. Actins are highly conserved proteins that are involved in cell motility, structure, and integrity. This actin is a major constituent of the contractile apparatus and one of the two nonmuscle cytoskeletal actins.
基因ID	60
基因名	ACTB
Swiss	P60709
别名	ACTB;BRWS1;PS1TP5BP1;beta actin;actin

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn