

Snail Rabbit pAb

货号: B25125

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IF: HeLa cells , Homo sapiens , Rattus norvegicus WB: HeLa cells , Homo sapiens , Mus musculus , Rattus norvegicus , Paralichthys olivaceus ICC: Homo sapiens IHC: Paralichthys olivaceus
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	29kDa
实测分子量	29KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,293T,HeLa,Mouse kidney
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 200-267 of human SNAI1 (NP_005976.2).
序列	RTHTGEKPFSCPHCSRAFADRSNLR AHLQTHSDVKKYQCQACARTFSRMSLLHKHQESGCSGCPR

靶点信息

研究背景	The Drosophila embryonic protein snail is a zinc finger transcriptional repressor which downregulates the expression of ectodermal genes within the mesoderm. The nuclear protein encoded by this gene is structurally similar to the Drosophila snail protein, and is also thought to be critical for mesoderm formation in the developing embryo. At least two variants of a similar processed pseudogene have been found on chromosome 2.
基因ID	6615
基因名	SNAI1
Swiss	O95863
别名	SNAI1;SLUGH2;SNA;SNAH;SNAIL;SNAIL1;dJ710H13.1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn