

HNRNPDL Rabbit pAb

货号: B25251

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:100 IF/ICC: 1:50 - 1:100
理论分子量	27kDa/33kDa/46kDa
实测分子量	38kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,SW620,HepG2,Mouse brain,Mouse testis,Rat brain
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 261-420 of human HNR NPDL (NP_112740.1).
序列	NIELPMDTKTNERRGFCFITYTDEEPVKLLLESRYHQIGSGKCEIKVAQPKEVYRQQQQQKGGRGAAAGGRGGTRGRG RGQGQWNWQGFNNYYDQGYGNYSAYGGDQNYSGYGGYDYGNYGNYGYGQGYADYSGQSQSTYGKASRGGGNH QNNYQPY

靶点信息

研究背景	This gene belongs to the subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). The hnRNPs are RNA binding proteins and they complex with heterogeneous nuclear RNA (hnRNA). These proteins are associated with pre-mRNAs in the nucleus and appear to influence pre-mRNA processing and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some seem to shuttle between the nucleus and the cytoplasm. The hnRNP proteins have distinct nucleic acid binding properties. The protein encoded by this gene has two RRM domains that bind to RNAs. Three alternatively spliced transcript variants have been described for this gene. One of the variants is probably not translated because the transcript is a candidate for nonsense-mediated mRNA decay. The protein isoforms encoded by this gene are similar to its family member HNRPD.
基因ID	9987
基因名	HNRNPDL
Swiss	O14979
别名	HNRNPDL;HNRNP;HNRPD;JKTBP;JKTBP2;LGMD1G;laAUF1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn