

NDUFA7 Rabbit pAb

货号: **B25286**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	12kDa
实测分子量	13kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,HepG2,MCF-7,Mouse kidney,Mouse brain,Mouse liver,Mouse heart,Rat liver
细胞定位	Matrix side,Mitochondrion inner membrane,Peripheral membrane protein
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-113 of human NDUFA 7 (NP_004992.2).
序列	MASATRLIQRLRNWASGHDLQGKLQLRYQEISKRTQPPPKLPVGPSHKLSNNYYCTRDGRRESVPPSIIMSSQKALVSGKP AESSAVAATEKKAVTPAPPIKRWELSSDQPYL

靶点信息

研究背景	This gene encodes a subunit of NADH:ubiquinone oxidoreductase (complex I), which is a multiprotein complex located in the inner mitochondrial membrane. Complex I functions in the transfer of electrons from NADH to the respiratory chain.
基因ID	4701

基因名	NDUFA7
Swiss	O95182
别名	NDUFA7;B14.5a;CI-B14.5a

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn