

NEDD4L Rabbit pAb

货号: B25311

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	96-111kDa
实测分子量	113kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BT-474,Mouse brain,Rat brain
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-975 of human NEDD4 L (Q96PU5).
序列	MATGLGEPVYGLSEDEGESRILRVKVVSGIDLAKKDIFGASDPYVKLSLYVADENRELALVQTKTIKKTLPKWNEEFYFRV NPSNHRLLFEVFDENRLTRDDFLGQVDVPLSHLPTEDPTMERPYTFKDFLLRPRSHKSRVKGFLRLKMAYMPKNGGQDEE NSDQRDDMEHGWEVVDSDNSASQHQEELPPPLPPGWEEKVDNLGRYYVNHNNRRTTQWHRPSLMDVSSSEDNNIR QINQEAHRRFRSRRHISEDLEPEPSEGGDVPEPWETISEEVNIAGDSLGLALPPPPASPGSRTSPQELSEELSRRLQITPDS NGEQFSSLIQREPSSRLRSCSVTDAVAEQGHLPPPSAPAGRARSSTVTGGEEPTPSVAYVHTTTPGLPSGWEERKDAKGR YYVNHNNRRTTTWTRPIMQLAEDGASGSATNSNNHLIEPQIRRPRLSSPTVTLSAPLEGAKDSPVRRRAVKDTLSNPQSPQP SPYNSPKPQHKTQSFLPPGWEMRIAPNGRPFFIDHNTKTTTWEDPRLKFPVHMRSKTSNPNLDGLPLPPGWEERIHLDG RTFYIDHNSKITQWEDPRLQNPAITGPAVPYSREFKQKYDYFRKKLKKPADIPNRFEMKLHRRNNIFEESYRRIMSVKRPDVLK ARLWIEFESEKGLDYGGVAREWFFLLSKEMFNPYYGLFEYSATDNYTLQINPNSGLCNEDHLSYFTFIGRVAGLAVFHGKLL DGFIRPFYKMMLGKQITLNDMESVDSEYYNSLKWILENDPTELDLMFIDEENFGQTYQVDLKPNGSEIMVTNENKREYI DLVIQWRFVNRVQKQMNAFLGFTTELLPIDLIKFDENELELLMCGLDGVDVNDWRQHSIYKNGYCPNHPVIQWFWKAVL LMDAEKRIRLLQFVTGTSRVPMNGFAELYGSNGPQLFTIEQWGSPEKLPRAHTCFNRDLPPYETFEDLREKLLMAVENAQ GFEGVD

靶点信息

研究背景	This gene encodes a member of the Nedd4 family of HECT domain E3 ubiquitin ligases. HECT domain E3 ubiquitin ligases transfer ubiquitin from E2 ubiquitin-conjugating enzymes to protein substrates, thus targeting specific proteins for lysosomal degradation. The encoded protein mediates the ubiquitination of multiple target substrates and plays a critical role in epithelial sodium transport by regulating the cell surface expression of the epithelial sodium channel, ENaC. Single nucleotide polymorphisms in this gene may be associated with essential hypertension. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
基因ID	23327
基因名	NEDD4L
Swiss	Q96PU5
别名	NEDD4L;NEDD4-2;NEDD4.2;PVNH7;RSP5;hNEDD4-2

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn