

TFE3 Rabbit pAb

货号: **B25315**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:100 - 1:500
理论分子量	11kDa/61kDa
实测分子量	75KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,U2OS,NIH/3T3,C6
细胞定位	Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human TFE3 (N P_006512.2).
序列	MSHAAEPARDGVEASAEGPRAVFVLLLEERRPADSAQLLSLNPESGIVADIELENVLDPDSFYELKSQPLPLRSSLPISLQ ATPATPATLSASSSAGGSRTPAMSSSSSSRVLLRQQLMRAQAQEQERRERREQAAAAPFPSPAPASPAISVVGVSAG

靶点信息

研究背景	This gene encodes a basic helix-loop-helix domain-containing transcription factor that binds MUE3-type E-box sequences in the promoter of genes. The encoded protein promotes the expression of genes downstr eam of transforming growth factor beta (TGF-beta) signaling. This gene may be involved in chromosomal translocations in renal cell carcinomas and other cancers, resulting in the production of fusion proteins. T ranslocation partners include PRCC (papillary renal cell carcinoma), NONO (non-POU domain containing, o ctamer-binding), and ASPSCR1 (alveolar soft part sarcoma chromosome region, candidate 1), among oth er genes. Alternative splicing results in multiple transcript variants.
------	--

基因ID	7030
基因名	TFE3
Swiss	P19532
别名	TFE3;RCCP2;RCCX1;TFEA;bHLHe33

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn