

WNT9A Rabbit pAb

货号: B25316

产品信息

| | |
|-------|--|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | |
| 应用 | WB |
| 推荐浓度 | WB: 1:500 - 1:2000 |
| 理论分子量 | 40kDa |
| 实测分子量 | 30kDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | Mouse lung |
| 细胞定位 | Secreted,extracellular matrix,extracellular space |
| 纯化 | Affinity purification |

抗原信息

| | |
|------|--|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 250-350 of human WNT 9A (NP_003386.1). |
| 序列 | ALKVGSTTNEAAGEAGAI SP PRGRASGAGGSDPLRTPELVHLDDSPSFCLAGRFSPGTAGRRCHREKNCESICCGRGHN TQSRVVTRPCQCQVRWCCYVE |

靶点信息

| | |
|------|---|
| 研究背景 | The WNT gene family consists of structurally related genes that encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. This gene is a member of the WNT gene family. It is expressed in gastric cancer cell lines. The protein encoded by this gene shows 75% amino acid identity to chicken Wnt14, which has been shown to play a central role in initiating synovial joint formation in the chick limb. This gene is clustered with another family member, WNT3A, in the chromosome 1q42 region. |
|------|---|

| | |
|-------|-------------|
| 基因ID | 7483 |
| 基因名 | WNT9A |
| Swiss | O14904 |
| 别名 | WNT9A;WNT14 |

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn