

# RNASEL Rabbit pAb

货号: **B25405**

产品信息

|       |                                                                                                        |
|-------|--------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                        |
| 宿主    | Rabbit                                                                                                 |
| 克隆性   | Polyclonal                                                                                             |
| 预测反应  | <b>WB:</b> Homo sapiens                                                                                |
| 应用    | <a href="#">WB</a>                                                                                     |
| 推荐浓度  | <b>WB:</b> 1:100 - 1:500                                                                               |
| 理论分子量 | 73kDa/83kDa                                                                                            |
| 实测分子量 | 80KDa                                                                                                  |
| 形式    | Liquid                                                                                                 |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                           |
| 阳性对照  | THP-1,Mouse spleen,Rat spleen                                                                          |
| 细胞定位  | Cytoplasm,Mitochondrion                                                                                |
| 纯化    | Affinity purification                                                                                  |

抗原信息

|      |                                                                                                                                                                                                                                                                                            |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 472-741 of human RNA SEL (NP_066956.1).                                                                                                                                                                      |
| 序列   | VQELHLSCGYTHQDLQPQNILIDSKKAAHLADFDKSIKWAGDPQEVKRDLEDLGRVLVLYVVKKGSISFEDLKAQSNEEVV<br>QLSPDEETKDLIHLRFHPGHEVRDCLSDLLGHPFFWTWESRYRTLNRVGNESDIKTRKSESEILRLQPGPSEHSKSFDKW<br>TTKINECVMKKMNKFYEKRGNFYQNTVGDLLKFIRNLGEHIDEEKHKMKLKIGDPSLYFQKTFPDLVIYVYTKLQNTTEYRK<br>HFPQTHSPNKPQCDGAGGASGLASPGC |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                       |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a component of the interferon-regulated 2-5A system that functions in the antiviral and antiproliferative roles of interferons. Mutations in this gene have been associated with predisposition to prostate cancer and this gene is a candidate for the hereditary prostate cancer 1 (HPC1) allele. |
| 基因ID | 6041                                                                                                                                                                                                                                                                                                                  |

|       |                   |
|-------|-------------------|
| 基因名   | RNASEL            |
| Swiss | Q05823            |
| 别名    | RNASEL;PRCA1;RNS4 |

产品验证

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)