

IRAK4 Rabbit pAb

货号: **B25427**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	37kDa/51kDa
实测分子量	55KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat lung
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human IRAK4 (NP_057207.2).
序列	MNKPITPSTYVRCLNVGLIRKLSDFIDPQEGWKKLAVAIAKKPSGDDRYNQFHIRRFEALLQTGKSPTSELLFDWGTTNCTV GDLVDLLIQNEFFAPASLLLPDAVPKTANTLPSKEAITVQKQMPFCDKDRTLMTVPQNLEQSYMPPDSSSPENKSLEVSD TRFHFSFSFYELKNVTNNFDERPISVGGNKMGEFGFVV

靶点信息

研究背景	This gene encodes a kinase that activates NF-kappaB in both the Toll-like receptor (TLR) and T-cell recept or (TCR) signaling pathways. The protein is essential for most innate immune responses. Mutations in this gene result in IRAK4 deficiency and recurrent invasive pneumococcal disease. Multiple transcript variants encoding different isoforms have been found for this gene.
------	---

基因ID	51135
基因名	IRAK4
Swiss	Q9NWZ3
别名	IRAK4;IPD1;IRAK-4;NY-REN-64;REN64

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn