

# Phospho-Jak3-Y980/981 Rabbit pAb

货号: **B25429**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>WB:</b> Rattus norvegicus
应用	<a href="#">WB</a>
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	68kDa/121kDa/125kDa
实测分子量	110kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	Cytoplasm,Endomembrane system,Peripheral membrane protein
纯化	Affinity purification

抗原信息

抗原信息	A synthetic phosphorylated peptide around Y980 & Y981 of human Jak3 (NP_000206.2).
序列	KDYVVV

靶点信息

研究背景	The protein encoded by this gene is a member of the Janus kinase (JAK) family of tyrosine kinases involved in cytokine receptor-mediated intracellular signal transduction. It is predominantly expressed in immune cells and transduces a signal in response to its activation via tyrosine phosphorylation by interleukin receptors. Mutations in this gene are associated with autosomal SCID (severe combined immunodeficiency disease).
基因ID	3718

基因名	JAK3
Swiss	P52333
别名	JAK3;JAK-3;JAK3_HUMAN;JAKL;L-JAK;LJAK

产品验证

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)