

TYK2 Rabbit pAb

货号: B25438

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus,Homo sapiens
应用	WB IHC IF/ICC
推荐浓度	WB: 1:100 - 1:500 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	133kDa
实测分子量	134KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Raji
细胞定位	cytoplasm,cytoskeleton,cytosol,extracellular exosome,nucleus
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1102-1182 of human TYK2 (NP_003322.3).
序列	QSPPTKFLELIGIAQGQMTVLRLTELLERGERLPRDKCPCEVYHLMKNCWETEASFRPTFENLIPILKTVHEKYQGQAPS

靶点信息

研究背景	This gene encodes a member of the tyrosine kinase and, more specifically, the Janus kinases (JAKs) protein families. This protein associates with the cytoplasmic domain of type I and type II cytokine receptors and promulgate cytokine signals by phosphorylating receptor subunits. It is also component of both the type I and type III interferon signaling pathways. As such, it may play a role in anti-viral immunity. A mutation in this gene has been associated with hyperimmunoglobulin E syndrome (HIES) - a primary immunodeficiency characterized by elevated serum immunoglobulin E.
------	---

基因ID	7297
基因名	TYK2
Swiss	P29597
别名	TYK2;IMD35;JTK1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn