

Phospho-IKKα/β-S176/180 Rabbit pAb

货号: **B25489**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens , Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	84kDa/29kDa/79kDa/85kDa/86kDa
实测分子量	85KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	cytoplasm,cytoplasmic side of plasma membrane,cytosol,IkappaB kinase complex,nucleoplasm
纯化	Affinity purification

抗原信息

抗原信息	A synthetic phosphorylated peptide around S176 & S180 of human IKK-alphaCHUK (NP_001269.3).
序列	SLCTS

靶点信息

研究背景	This gene encodes a member of the serine/threonine protein kinase family. The encoded protein, a component of a cytokine-activated protein complex that is an inhibitor of the essential transcription factor NF-kappa-B complex, phosphorylates sites that trigger the degradation of the inhibitor via the ubiquitination pathway, thereby activating the transcription factor. [provided by RefSeq, Jul 2008]
基因ID	1147,3551
基因名	CHUK,IKBKB

Swiss	O15111,O14920
别名	

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn