

PEG400

货号: **B26660**

PEG400

产品信息

生物活性	PEG400 is a strongly hydrophilic polyethylene glycol used as an excellent solvent for a large number of substances. PEG400 is widely used in a variety of pharmaceutical formulations.
CAS	
中文名称	
分子量	
体外研究	
体内研究	
形式	Liquid
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	
溶解性	In Vitro: H₂O : 100 mg/mL (Need ultrasonic) DMSO : 100 mg/mL (Need ultrasonic)
纯度	