

Alkaline phosphatase

货号: **B26692**

Alkaline phosphatase

产品信息

生物活性	Alkaline phosphatase (Apase) is a membrane-bound glycoprotein that catalyzes the hydrolysis of phosphate monoesters at basic pH values. Alkaline phosphatase can be used for molecular biology and enzyme immunoassay.
CAS	9001-78-9
中文名称	碱性磷酸酶
分子量	
体外研究	
体内研究	
形式	Solid
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	In Vitro: H₂O : 100 mg/mL (Need ultrasonic) In Vivo: 请根据您的实验动物和给药方式选择适当的溶解方案。以下溶解方案都请先按照 In Vitro 方式配制澄清的储备液，再依次添加助溶剂： ——为保证实验结果的可靠性，澄清的储备液可以根据储存条件，适当保存；体内实验的工作液，建议您现用现配，当天使用； 以下溶剂前显示的百分比是指该溶剂在您配制终溶液中的体积占比；如在配制过程中出现沉淀、析出现象，可以通过加热和/或超声的方式助溶 1. 请依序添加每种溶剂： PBS Solubility: 50 mg/mL (Infinity mM); Clear solution; Need ultrasonic *以上所有助溶剂都可在 MCE 网站选购。
纯度	