

Nattokinase

货号: **B26698**

Nattokinase

产品信息

生物活性	Nattokinase is a potent fibrinolytic enzyme. Nattokinase can break down blood clots by directly hydrolyzing fibrin and plasmin substrate. Nattokinase can be used for the research of cardiovascular diseases.
CAS	133876-92-3
中文名称	
分子量	
体外研究	<p>Nattokinase dissolves 94% blood clot at ten minutes and degrades the fibrin clot after two hours of incubation.</p> <p>The accuracy of these methods have not been independently confirmed. They are for reference only.</p>
体内研究	
形式	Solid
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	<p>In Vitro: H₂O : 100 mg/mL (Need ultrasonic)</p> <p>In Vivo: 请根据您的实验动物和给药方式选择适当的溶解方案。以下溶解方案都请先按照 In Vitro 方式配制澄清的储备液，再依次添加助溶剂：</p> <p>——为保证实验结果的可靠性，澄清的储备液可以根据储存条件，适当保存；体内实验的工作液，建议您现用现配，当天使用； 以下溶剂前显示的百分比是指该溶剂在您配制终溶液中的体积占比；如在配制过程中出现沉淀、析出现象，可以通过加热和/或超声的方式助溶</p> <ul style="list-style-type: none">1. 请依序添加每种溶剂： PBS <p>Solubility: 50 mg/mL (Infinity mM); Clear solution; Need ultrasonic</p> <p>*以上所有助溶剂都可在 MCE 网站选购。</p>
纯度	

