


Polyvinylpyrrolidone

货号: **B26700**

产品信息

生物活性	Polyvinylpyrrolidone is a compound which has been widely tested and used in human and veterinary medicine as an effective wound healing accelerator and disinfectant when combined with iodine and other compounds.
CAS	9003-39-8
中文名称	聚乙烯吡咯烷酮
分子量	111.1418 (monomer)
体外研究	
体内研究	
形式	Solid
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	
溶解性	In Vitro: H₂O : ≥ 50 mg/mL DMSO : 25 mg/mL (Need ultrasonic) *"≥" means soluble, but saturation unknown.
纯度	≥98.0%