


# 4-Aminobenzaldehyde

货号: **B26726**

## 产品信息

生物活性	4-Aminobenzaldehyde (p-aminobenzaldehyde) is a useful synthetic reagent and monomer that can be used to synthesize monoazo dyes and photocurable ion exchange resins. 4-Aminobenzaldehyde is also a corrosion inhibitor of metals.
CAS	556-18-3
中文名称	对氨基苯甲醛
分子量	121.14
体外研究	
体内研究	
形式	Solid
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	4°C, protect from light, stored under nitrogen
溶解性	In Vitro: <b>DMSO : (insoluble or slightly soluble)</b>
纯度	≥98.0%