


4-Dodecylbenzenesulfonic acid

货号: **B26832**



产品信息

生物活性	4-Dodecylbenzenesulfonic acid is a biochemical reagent that can be used as a biological material or organic compound for life science related research.																
CAS	121-65-3																
中文名称	十二烷基苯磺酸																
分子量	326.49																
体外研究																	
体内研究																	
形式	Liquid																
运输条件	Room temperature in continental US; may vary elsewhere.																
保存条件																	
溶解性	In Vitro: DMSO : ≥ 100 mg/mL (306.29 mM) *"≥" means soluble, but saturation unknown. 配制储备液 <table><tr><th>浓度</th><th>溶剂</th><th>体积</th><th>质量</th></tr><tr><td>1 mM</td><td>3.0629 mL</td><td>15.3144 mL</td><td>30.6288 mL</td></tr><tr><td>5 mM</td><td>0.6126 mL</td><td>3.0629 mL</td><td>6.1258 mL</td></tr><tr><td>10 mM</td><td>0.3063 mL</td><td>1.5314 mL</td><td>3.0629 mL</td></tr></table> * 请根据产品在不同溶剂中的溶解度选择合适的溶剂配制储备液；一旦配成溶液，请分装保存，避免反复冻融造成的产品失效。 储备液的保存方式和期限：-80°C, 6 months; -20°C, 1 month。-80°C 储存时，请在 6 个月内使用，-20°C 储存时，请在 1 个月内使用。	浓度	溶剂	体积	质量	1 mM	3.0629 mL	15.3144 mL	30.6288 mL	5 mM	0.6126 mL	3.0629 mL	6.1258 mL	10 mM	0.3063 mL	1.5314 mL	3.0629 mL
浓度	溶剂	体积	质量														
1 mM	3.0629 mL	15.3144 mL	30.6288 mL														
5 mM	0.6126 mL	3.0629 mL	6.1258 mL														
10 mM	0.3063 mL	1.5314 mL	3.0629 mL														
纯度																	