


Pseudouridine 5'-triphosphate

货号: B26887



产品信息

生物活性	pseudouridine-5'-triphosphate (Pseudo-UTP) is one of the most commonly used modified nucleoside for the polymerase-mediated synthesis of RNA molecules. Compared with uridine-containing unmodified mRNAs, the application of pseudouridine-containing modified mRNAs exhibits better nuclease stability, immunogenicity, and translational properties..																		
CAS	1175-34-4																		
中文名称																			
分子量	484.14																		
体外研究																			
体内研究																			
形式																			
运输条件	Room temperature in continental US; may vary elsewhere.																		
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.																		
溶解性	<p>In Vitro:</p> <p>H₂O : 250 mg/mL (516.38 mM; Need ultrasonic)</p> <p>配制储备液</p> <table><tr><th>浓度</th><th>溶剂</th><th>体积</th><th>质量</th></tr><tr><td>1 mg</td><td></td><td></td><td></td></tr><tr><td>5 mg</td><td></td><td></td><td></td></tr><tr><td>10 mg</td><td></td><td></td><td></td></tr></table> <p>1 mM 2.0655 mL10.3276 mL20.6552 mL</p> <p>5 mM 0.4131 mL2.0655 mL 4.1310 mL</p> <p>10 mM0.2066 mL1.0328 mL 2.0655 mL</p>			浓度	溶剂	体积	质量	1 mg				5 mg				10 mg			
浓度	溶剂	体积	质量																
1 mg																			
5 mg																			
10 mg																			
纯度																			