


# 5,10,15,20-Tetrakis(p-tolyl)porphyrin

货号: B26897



产品信息

|      |                                                                                                                                                                 |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 生物活性 | 5,10,15,20-Tetrakis(p-tolyl)porphyrin is a biochemical reagent that can be used as a biological material or organic compound for life science related research. |
| CAS  | 14527-51-6                                                                                                                                                      |
| 中文名称 | 5 10 15 20-四对甲苯基-21H 23H-卟吩                                                                                                                                     |
| 分子量  | 670.84                                                                                                                                                          |
| 体外研究 |                                                                                                                                                                 |
| 体内研究 |                                                                                                                                                                 |
| 形式   | Solid                                                                                                                                                           |
| 运输条件 | Room temperature in continental US; may vary elsewhere.                                                                                                         |
| 保存条件 | 4°C, protect from light, stored under nitrogen                                                                                                                  |
| 溶解性  |                                                                                                                                                                 |
| 纯度   |                                                                                                                                                                 |