


Biotin-C1-PEG3-C3-amido-C5-Gly-Arg-Gly-N3 TFA

货号: B26946



产品信息

生物活性	Biotin-C1-PEG3-C3-amido-C5-Gly-Arg-Gly-N3 TFA is used for detection of modification site for N-myristoylated and GPI-anchored proteins in blood-stage <i>P. falciparum</i> .																
CAS																	
中文名称																	
分子量	970.07																
体外研究																	
体内研究																	
形式	Solid																
运输条件	Room temperature in continental US; may vary elsewhere.																
保存条件	-20°C, sealed storage, away from moisture																
溶解性	In Vitro: DMSO : 100 mg/mL (103.09 mM; Need ultrasonic) 配制储备液 <table><tr><td>浓度</td><td>溶剂</td><td>体积</td><td>质量</td></tr><tr><td>1 mg</td><td></td><td></td><td></td></tr><tr><td>5 mg</td><td></td><td></td><td></td></tr><tr><td>10 mg</td><td></td><td></td><td></td></tr></table> 1 mM 1.0309 mL5.1543 mL10.3085 mL 5 mM 0.2062 mL1.0309 mL2.0617 mL 10 mM0.1031 mL0.5154 mL1.0309 mL * 请根据产品在不同溶剂中的溶解度选择合适的溶剂配制储备液；一旦配成溶液，请分装保存，避免反复冻融造成的产品失效。 储备液的保存方式和期限：-80°C, 6 months; -20°C, 1 month (sealed storage, away from moisture)。-80°C 储存时，请在 6 个月内使用，-20°C 储存时，请在 1 个月内使用。	浓度	溶剂	体积	质量	1 mg				5 mg				10 mg			
浓度	溶剂	体积	质量														
1 mg																	
5 mg																	
10 mg																	
纯度	≥98.0%																