


# (E)-Hex-3-enoic acid

货号: **B26948**

## 产品信息

|      |                                                                                                                                                |
|------|------------------------------------------------------------------------------------------------------------------------------------------------|
| 生物活性 | (E)-Hex-3-enoic acid is a biochemical reagent that can be used as a biological material or organic compound for life science related research. |
| CAS  | 1577-18-0                                                                                                                                      |
| 中文名称 | 反-3-己烯酸                                                                                                                                        |
| 分子量  | 114.14                                                                                                                                         |
| 体外研究 |                                                                                                                                                |
| 体内研究 |                                                                                                                                                |
| 形式   | 12°C Liquid                                                                                                                                    |
| 运输条件 | Room temperature in continental US; may vary elsewhere.                                                                                        |
| 保存条件 |                                                                                                                                                |
| 溶解性  |                                                                                                                                                |
| 纯度   |                                                                                                                                                |