


Phosphoric acid (sodium hydrate), ≥99.0%

货号: B26997



产品信息

生物活性	Phosphoric acid (sodium hydrate), ≥99.0% is used as a buffer capacity reagent in molecular biology , biochemistry, and chromatography. For the preparation of biological buffers. It is also used to purify antibodies , as a laxative, and in combination with other sodium phosphates
CAS	13472-35-0
中文名称	磷酸二氢钠(二水)
分子量	156.01
体外研究	
体内研究	
形式	Solid
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	4°C, sealed storage, away from moisture
溶解性	
纯度	