

alpha-(2,3)-Sialyltransferase(PmST1)

货号: B27005

alpha-(2,3)-Sialyltransferase(PmST1)

产品信息

生物活性	alpha-(2,3)-Sialyltransferase (PmST1), a sialyltransferase, can be isolated from <i>Pasteurella multocida</i> . alpha-(2,3)-Sialyltransferase can be converted into regioselective α 2-6-sialyltransferase by saturation mutagenesis and regioselective screening.
CAS	71124-51-1
中文名称	
分子量	
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	