

Ficain

货号: B27016

Ficain

产品信息

生物活性	Ficain is an enzyme extract composed of several proteases that can be isolated from <i>Ficus hispida</i> L. and the latex of fig (<i>Ficus carica</i>). Ficain has different specificities in different proportions during fruit ripening . Ficain is widely used in protein hydrolysis, food, production of bioactive peptides and antibody fragments.
CAS	9001-33-6
中文名称	
分子量	
体外研究	<p>Ficain (0-8 µg/mL) 以酪蛋白为底物，在 280 nm 处测量三氯乙酸滤液的吸光度，检测 Ficain 活性。</p> <p>Ficain (0-200 µg/mL) 在 p-nitrophenyl carbobenzoxyglycinate 为底物时，以剂量依赖性方式增加 400 nm 处的吸光度。</p> <p>The accuracy of these methods have not been independently confirmed. They are for reference only.</p>
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	