


G3-CNP

货号: B27081



产品信息

生物活性	G3-CNP is an α -amylase substrate. The absorbance of G3-CNP cleavage product 2-chloro-4-nitrophenol is measured at 405 nm, which can be used to detect enzyme activity.
CAS	118291-90-0
中文名称	
分子量	659.98
体外研究	G3-CNP (50 μ L; 10 min) 在 37°C 条件下, 被 α -淀粉酶裂解产生 2-氯-4-硝基苯酚。 The accuracy of these methods have not been independently confirmed. They are for reference only.
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	