

# Hyaluronan synthase

货号: B27158

## Hyaluronan synthase

### 产品信息

生物活性	Hyaluronan synthase is a membrane protein that requires only Mg <sup>+2</sup> and two sugar-UDP substrates (GlcU A-UDP and GlcNAc-UDP) to polymerize HA chains. Hyaluronan synthase catalyzes the biosynthesis of hyal uronic acid (HA).
CAS	39346-43-5
中文名称	
分子量	
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	