


# Tetramethylammoniump-toluenesulfonate

货号: **B27349**

## 产品信息

|      |                                                                                                                                                                 |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 生物活性 | Tetramethylammoniump-toluenesulfonate is a biochemical reagent that can be used as a biological material or organic compound for life science related research. |
| CAS  | 3983-91-3                                                                                                                                                       |
| 中文名称 | 四甲基对甲苯磺酸铵 25G                                                                                                                                                   |
| 分子量  | 245.34                                                                                                                                                          |
| 体外研究 |                                                                                                                                                                 |
| 体内研究 |                                                                                                                                                                 |
| 形式   | Solid                                                                                                                                                           |
| 运输条件 | Room temperature in continental US; may vary elsewhere.                                                                                                         |
| 保存条件 |                                                                                                                                                                 |
| 溶解性  |                                                                                                                                                                 |
| 纯度   |                                                                                                                                                                 |