

D-Lactate dehydrogenase

货号: B27355

D-Lactate dehydrogenase

产品信息

生物活性	D-Lactate dehydrogenase (D-LDH) is an oxidoreductase that uses NAD ⁺ or NADP ⁺ as an acceptor and acts on the donor CH-OH group, and can catalyze the oxidation of D-lactate to pyruvate. D-Lactate dehydrogenase widely exists in bacteria and fungi, and is often used in biochemical research.
CAS	9028-36-8
中文名称	
分子量	
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	