

Triosephosphate isomerase

货号: **B27421**

Triosephosphate isomerase

产品信息

生物活性	Triosephosphate isomerase (TPI) is a glycolytic enzyme. Triosephosphate isomerase fast interconverts dihydroxyacetone phosphate and D: -glyceraldehyde-3-phosphate, its catalytic site is at the dimer interface.
CAS	9023-78-3
中文名称	
分子量	
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	